


EDT Syringe

Enhanced Dosing Technology

REF: ZRS02EG-CS



SJJ Solutions' Enhanced Dosing Technology (EDT) syringe revolutionizes precise dosing by enabling accurate measurement of volumes smaller than 50µL. By using the red dosing line instead of the piston, the syringe enables precise dosing in the 10µL to 200µL range.

Together with its low residual volume and optimal return, the EDT syringe offers accurate and smooth dosage control, complemented by a favourable volumetric flow rate for consistent and reliable performance.

The syringe is easy-to-use and offers great advantages for injecting sub-50µL.

Specifications

- Low residual volume & optimal accuracy
- Provides a smooth, even glide throughout the injection
- Manufactured without silicone oil and is USP 789 compliant
- User friendly & ergonomic design
- Favorable volumetric flow rate
- Anti-VEGF stability data (up to 30 days)
- Luer lock connection
- CE & FDA approved

Research

[Accuracy and residual volume of intravitreal injection syringes in administration of very small volumes](#)